

RH
scientific

Tissue Flotation Bath
TFB-1207



RH SCIENTIFIC

Meet the RH Scientific

TISSUE FLOTATION BATH



- **TEFLON COATED TUB**
- **USER FRIENDLY**
- **MICROPROCESSOR CONTROL**
- **DIGITAL DISPLAY**
- **SILICON HEATER**
- **CONSTANT TEMPERATURE**

for any query contact :

SINGH ENTERPRISES

FARIDABAD ,HARYANA ,121001

Phone : 9319193775,8448331034

Email- singhenterprises469@gmail.com

www.rhscientific.com

PRODUCT DESCRIPTION:

TISSUE FLOTATION BATH TFB-1207 is a simple, economical, designed to lay section of Histological samples. It provides high thermal conductivity rates and outstanding scratches resistance due to its special coating overheating protection system along with a stand- by mode.

TECHNICAL SPECIFICATION:

- ❖ Precise temperature control.
- ❖ Built in memory for storing set temperatures.
- ❖ Jet-Black surface for easy identification of Tissues.
- ❖ Suitable to work on
220V, Single phase,
50Hz AC supply.
 - Water Temperature Range: Ambient to 99°C
 - Temperature Tolerance: +1°C
 - Power: 200V, 50Hz
 - Tub Dimensions: 300x200x60mm

for any query contact :

SINGH ENTERPRISES

FARIDABAD ,HARYANA ,121001

Phone : 9319193775,8448331034

Email- singhenterprises469@gmail.com