

# Tissue Block Cabinet for 20000 Block



Tissue Block/Casette/Ring Storage Cabinet and drawers made up of powder coated mildsteel.

Each drawer is having compartments to store the tissue cassette/ring/block (back to back in close pack manner).

Cardplate and Puller knob fixed on front of the drawer. Cabinet provided with hinged door with locking system and handles.

Available capacity: 20000

## Dimensions:

Cabinet Unit: 36 inch (H) x 17.75inch(D) x 60 (w)

Drawer size: 18.5x 14.5inch

*Brass pull knobs and card holders on each dual slotted drawer. Included foam inserts keep cassette neatly standing.*

*Interlocking system makes it possible to safely store up to 20,000 cassette on a single stack.*



For any query contact us:

**Singh enterprises**

Faridabad ,Haryana

9319134775,8448331034

Singhenterprises469@gmail.com