

RH scientific

HOT PLATE TSW-0709



RH SCIENTIFIC

Meet the RH Scientific

HOT PLATE TSW-0709



- **CONSTANT TEMPERATURE**
- **TEFLON COATED TOP**
- **USER FRIENDLY MICROPROCESSOR CONTROL**
- **SILICON HEATER**
- **DIGITAL DISPLAY**

for any query contact :

SINGH ENTERPRISES

FARIDABAD ,HARYANA ,121001

Phone : 9319193775,8448331034

Email- singhenterprises469@gmail.com

www.rhscientifics.com

PRODUCT DESCRIPTION:

HOT PLATE is used in histopathology laboratory for drying slides at a uniform temperature just below the melting point of wax so as to provide flattening without damage to the section. This HOT PLATE is fabricated from Top notch materials as pre- defined industry Standards. Our Slide Warming Tables are known for their optimum utility and high performance.

TECHNICAL SPECIFICATION:

- ❖ Feather touch keypad with LCD display showing:
Date, Time, Day, Set value, Present value.
- ❖ Microprocessor Precise Temperature controls.
Built in memory for storing temperature.
- ❖ Uniform temperature throughout the surface.
- ❖ Corrosion Resistant material for durability.
 - Temperature Range: ambient to 100°C
 - Temperature Tolerance: +/-1°C
 - Slide Warmer Dimensions: 340x240mm
 - Power: 220V 50Hz
 - Overall dimensions: 360x360x210mm

for any query contact :

SINGH ENTERPRISES

FARIDABAD ,HARYANA ,121001

Phone : 9319193775,8448331034

Email- singhenterprises469@gmail.com